



Bi-level PAP With LCD Screen

Specification

LCD display

With pressure and flow waveform display

With real-time monitoring 8 kinds of respiratory parameters

With independent medical grade humidifier

Trolley with patch board, oxygen cylinder slot

With full face mask and SD card

Mode: CPAP, S, T, ST, VGPS

Pressure: 4-30 cm H₂O

Trigger adjustment

Cycle adjustment

Slope adjustment

Noise: ≤ 32 dB(A)

Certification: CE /CFDA approved

Package Size : 1 CTN 63×63×105cm

N.W.(kg): 25

G.W.(kg): 30







